
SEAL STAIN ABLE



A MORE SUSTAINABLE PACKAGING IN TWO STEPS

hawa







Step 01 | Sustainable SBS Sealers



GreenTek sealers by hawo such as the hm 950 DC-VI NanoPak are equipped with sustainable and resource-saving technologies. These devices only need a fraction of the energy required by conventional systems, the carbon-footprint is correspondingly low.

hawo sealing devices are of course produced with self-generated solar power and shipped in plasticfree transport packaging, completely in line with our ConceptZero.







Step 02 | Plastic-free Protective Packaging



ORGANIXprotect is the first **sustainable protective packaging** for sterile barrier systems. As a so-called Dust Cover, it fulfills the requirements for protective packaging according to DIN EN ISO 11607-01 by providing a physical protection for the SBS and its contents.

ORGANIXprotect is made exclusively from sustainable raw materials and is **100% plastic-free, compostable and recyclable.**





hm 950 DC-VI Nanopak

A state-of-the-art heating system brings NanoPak and other models to their operating temperature in a record time and operates with a highly efficient 24 volt system that reduces power consumption to a minimum.



hd 480 WSI-V ValiPakTOUCH

The stand-by function or automatic switch-off when not in use is a standard feature of more and more hawo devices. It perfectly completes the GreenTek concept without any compromises regarding the quality of the sealing results.







Rolls:

- 200 mm x 400 m
- 250 mm x 400 m
- 300 mm x 400 m

ORGANIX protect R

Standard rolls in various widths for an individual adaptation to the size of each sterile barrier system. With the help of a film sealing device, the protective packaging can be produced precisely and cost-effectively on site.



Self-adhesive pouches:

- 200 x 300 + 50 mm
- 250 x 400 + 50 mm
- 300 x 500 + 50 mm

Tolerance range +/- 2 mm. Other sizes on request.

ORGANIX protect S

Prefabricated self-adhesive pouches in different sizes for an even more convenient creation of the protective packaging – ideally suited for CSSDs with limited space conditions or standardized SBS sizes.



One of our natural resources.

hawo GmbH Obere Au 2-4 74847 Obrigheim Germany

T +49 (0) 6261/9770-0 F +49 (0) 6261/9770-69

info@hawo.com www.hawo.com

Important advice:

Non-sterile. Only pack dry and cool products/ sterile goods, risk of condensation.

ORGANIXprotect is not





ORGANIXprotect is not suitable as a safety bag for disinfected medical devices.